


# ORTHORA 200

EQUIPMENT

2506


CONSULTANT


DISTRIBUTOR

Operation	Voltage	Socket system	Color	Upholstery	Aspirator
<input type="checkbox"/> Right-handed person <input type="checkbox"/> Left-handed person	<input type="checkbox"/> 230V/50Hz <input type="checkbox"/> 100V/50Hz <input type="checkbox"/> 120V/60Hz	<input type="checkbox"/> Schuko <input type="checkbox"/> CH T13 <input type="checkbox"/> US/CAN <input type="checkbox"/> F/B	<input type="checkbox"/> RAL (Standard) RAL-No. _____  <input type="checkbox"/> Metallic (Surcharge) Description MIK _____	Upholstery Mikrona "T" <input type="checkbox"/> Standard <input type="checkbox"/> Memory Color: _____	<input type="checkbox"/> Dry aspirator syst. <input type="checkbox"/> Wet aspirator syst.  <input type="checkbox"/> None

**Re.Formance Line Working chair**



**Re.Forma Hi.**  
Amount: \_\_\_\_\_



**Re.Forma**  
Amount: \_\_\_\_\_




**Re.Forma Lite.**  
Amount: \_\_\_\_\_

**Standard**  
Black base


**Options**  
 Chromed base  
 Foot plate control for height adjustment (chromed)

**Upholstery**  
 Leather upholstery like unit  
 Leather upholstery \_\_\_\_\_



**Armrest**  
 Yes  
 No

**Attention!**  
 For safety reasons the armrest may not be used in combination with ORTHODESK 320.6545.xx .



**Instrument bracket**  
 (right or left mounted)  
 Yes  
 No

**Instrument bracket with**  
 Micromotor MX & MX2  
 Micromotor MC3 or turbine device  
 Ultrasonic EMS  
 Ultrasonic PIEZON  
 3-Function syringe LUZZANI  
 6-Function syringe LUZZANI



**ORTHORACK**

**ORTHORACK Back-of-head instrument (BOH)**

- BOH connected with OT200 unit
  - Air/water supply from OT200 unit
  - Air/water supply from BOH

**Autonomous water supply**



- Yes
- No



**Rack with aspirator**  
 Yes  OP-suction  
 Saliva suction  
 No

**Instrument rack (max. 4 instr.)**

- Micromotor MX2 LED
  - Instr. hose with water regulation
- Micromotor MC3 IR
- Micromotor MC3 LED
- Turbine device 4VLM LUX
- Ultrasonic PIEZON, EMS
- Ultrasonic NEWTRON, Satelec
- 3-Function syringe LUZZANI
- 6-Function syringe LUZZANI
- Polymerization lamp MINI LED ORTHO



**iPad-holder**

on floating table


- iPad mini 1 / 2 / 3 / 4
- iPad Air2 & Pro9,7 & 9,7"
- iPad Pro 10,5" & Pro 11" (2018) & 10,2"

**Multimedia**

- Yes
- No

**Adaptation**

- LCD-Screen 22" on mounting rod
- Intraoral camera



**OP lamp**

- Yes
- No

**Lamp**

- OP lamp LED MIKRONA (lamp case, RAL 9006)
  - sensor  switch
- OP lamp arm extension
- Ceiling OP lamp LED MIKRONA (Polaris)
  - RAL 9002  RAL 9010  RAL 9016
  - customer specific finish
  - sensor  switch
- ORTHOLUX 200 LED (ceiling model, color RAL 9002)



**Floating table**

- Yes
- No

**Floating table**

- without instrument bracket
- with instrument bracket VARIO

**Instrument bracket (max. 3 instruments)**


- Micromotor MX2 LED
  - Instr. hose with water regulation
- Micromotor MC3 IR
- Micromotor MC3 LED
- Turbine device 4VLM LUX
- Ultrasonic PIEZON, EMS
- Ultrasonic NEWTRON, Satelec
- 3-Function syringe LUZZANI
- 6-Function syringe LUZZANI
- Rack for micro blaster
  - Microetcher  Dento-Prep

**Equipment tray**

- Yes
- No


**Media supply**

- Media supply 1
- Media supply 2
- Media supply 3
- Water from bottle system




**Tray on comfort rack**

- Yes
- No



**Spittoon**

- with glass filler
  - Yes  Glass bowl (standard)
  - Ceramic bowl (optional)
- No



**Instruments positioned on chair base**

- Yes
- No

**Instruments**


- Micromotor MX2 LED
- Turbine device 4VLM LUX
- Air pistol FREUDING

**Quick-action coupling for external equipment**

- Air/water
  - Water from bottle system
- Air


**Water separation system**

- Autonomous water supply
  - 0.8lt Bottle
  - 1.2lt Bottle
- Metasys water aseptic & separation system WEK (apparatus beside the chair)



**Instrument flushing**

- CleanHub



**Rack with aspirator**

- Yes
  - OP-suction
  - Saliva suction
- No

**Instrument rack (max. 4 instruments)**

- Micromotor MX2 LED
  - Instr. hose with water regulation
- Micromotor MC3 IR
- Micromotor MC3 LED
- Turbine device 4VLM LUX
- Ultrasonic PIEZON, EMS
- Ultrasonic NEWTRON, Satelec
- 3-Function syringe LUZZANI
- 6-Function syringe LUZZANI
- Rack for micro blaster
  - with air connection  Microetcher
  - without air connection  Dento-Prep

**Instrument operation**


- 4-Function pedal DYN
- 5-Function pedal Pedal
  - built in (standard)
  - free standing (wireless)
  - free standing (cable connected)

**3-Function pedal for operating special functions**


- Yes
- No

**Operating functions**

- ORTHOLUX 200
- OP lamp ON/OFF
- special functions



**ORTHODESK**



**Media supply ORTHODESK**

- Yes (only possible, if no instrument is mounted on chair base)
- No

**Media supply (2 instruments)**

- Micromotor MX & MX2
- Micromotor MC3 IR
- Micromotor MC3 LED
- Turbine device 4VLM LUX
- 3-Function syringe LUZZANI
- 6-Function syringe LUZZANI

max. 1 instr.

Date:

Signature: