


ORTHORA 200

EQUIPMENT 2412

CONSULTANT


DISTRIBUTOR




ORTHORA 200

Operation	Voltage	Socket system	Color	Upholstery	Aspirator
<input type="checkbox"/> Right-handed person <input type="checkbox"/> Left-handed person	<input type="checkbox"/> 230V/50Hz <input type="checkbox"/> 100V/50Hz <input type="checkbox"/> 120V/60Hz	<input type="checkbox"/> Schuko <input type="checkbox"/> CH T13 <input type="checkbox"/> US/CAN <input type="checkbox"/> F/B	<input type="checkbox"/> RAL (Standard) RAL-No. _____ <input type="checkbox"/> Metallic (Surcharge) Description MIK _____	Upholstery Mikrona "T" <input type="checkbox"/> Standard <input type="checkbox"/> Memory Color: _____	<input type="checkbox"/> Dry aspirator syst. <input type="checkbox"/> Wet aspirator syst. <input type="checkbox"/> None

Re.Formance Line Working chair



☐ **Re.Forma Hi.**
Amount: _____



☐ **Re.Forma**
Amount: _____



☐ **Re.Forma Lite.**
Amount: _____


Standard
Black base

Options

- ☐ Chromed base
☐ Foot plate control for height adjustment (chromed)

Upholstery

- ☐ Leather upholstery like unit
☐ Leather upholstery _____




Armrest

- ☐ Yes
☐ No

Attention!

For safety reasons the armrest may not be used in combination with ORTHODESK 320.6545.xx .



Instrument bracket (right or left mounted)

- ☐ Yes
☐ No

Instrument bracket with

- ☐ Micromotor MX & MX2
☐ Micromotor MC3
 or turbine device
☐ Ultrasonic EMS
☐ Ultrasonic PIEZON
☐ 3-Function syringe LUZZANI
☐ 6-Function syringe LUZZANI



ORTHORACK

ORTHORACK Back-of-head instrument (BOH)

- ☐ BOH connected with OT200 unit
☐ Air/water supply from OT200 unit
☐ Air/water supply from BOH


Autonomous water supply

- ☐ Yes
☐ No




BOH

- ☐ standard
 (470mm deep)



BOH

- ☐ countersunk
 (570mm deep)



Rack with aspirator

- ☐ Yes ☒ OP-suction
 ☒ Saliva suction
☐ No

Instrument rack (max. 4 instr.)

- ☐ Micromotor MX2 LED
 ☐ Instr. hose with water regulation
☐ Micromotor MC3 IR
☐ Micromotor MC3 LED
☐ Turbine device 4VLM LUX
☐ Ultrasonic PIEZON, EMS
☐ Ultrasonic NEWTRON, Satelec
☐ 3-Function syringe LUZZANI
☐ 6-Function syringe LUZZANI
☐ Polymerization lamp MINI LED ORTHO

iPad-holder

on floating table


- ☐ iPad mini 1 / 2 / 3 / 4
- ☐ iPad Air2 & Pro9,7 & 9,7"
- ☐ iPad Pro 10,5" & Pro 11" (2018) & 10,2"

Multimedia

- ☐ Yes
- ☐ No

Adaptation

- ☐ LCD-Screen 22" on mounting rod
- ☐ Intraoral camera




OP lamp

- ☐ Yes
- ☐ No

Lamp

- ☐ OP lamp LED MIKRONA (lamp case, RAL 9006)
 - ☐ sensor ☐ switch
- ☐ OP lamp arm extension
- ☐ Ceiling OP lamp LED MIKRONA (Polaris)
 - ☐ RAL 9002 ☐ RAL 9010 ☐ RAL 9016
 - ☐ customer specific finish
 - ☐ sensor ☐ switch
- ☐ ORTHOLUX 200 LED (ceiling model, color RAL 9002)



Floating table


- ☐ Yes
- ☐ No

Floating table

- ☐ without instrument bracket
- ☐ with instrument bracket VARIO

Instrument bracket (max. 3 instruments)

- ☐ Micromotor MX2 LED
 - ☐ Instr. hose with water regulation
- ☐ Micromotor MC3 IR
- ☐ Micromotor MC3 LED
- ☐ Turbine device 4VLM LUX
- ☐ Ultrasonic PIEZON, EMS
- ☐ Ultrasonic NEWTRON, Satelec
- ☐ 3-Function syringe LUZZANI
- ☐ 6-Function syringe LUZZANI
- ☐ Rack for micro blaster
 - ☐ Microetcher ☐ Dento-Prep



Equipment tray


- ☐ Yes
- ☐ No

Media supply

- ☐ Media supply 1
- ☐ Media supply 2
- ☐ Media supply 3
- ☐ Water from bottle system


Tray on comfort rack

- ☐ Yes
- ☐ No




Water separation system

- ☐ Autonomous water supply
 - ☐ 0.8lt Bottle
 - ☐ 1.2lt Bottle
- ☐ Metasys water aseptic & separation system WEK (apparatus beside the chair)



Instrument flushing

- ☐ CleanHub



Spittoon


- with glass filler
 - ☐ Yes ☐ Glass bowl (standard)
 - ☐ Ceramic bowl (optional)
- ☐ No

Rack with aspirator

- ☐ Yes
 - ☒ OP-suction
 - ☒ Saliva suction
- ☐ No

Instrument rack (max. 4 instruments)

- ☐ Micromotor MX2 LED
 - ☐ Instr. hose with water regulation
- ☐ Micromotor MC3 IR
- ☐ Micromotor MC3 LED
- ☐ Turbine device 4VLM LUX
- ☐ Ultrasonic PIEZON, EMS
- ☐ Ultrasonic NEWTRON, Satelec
- ☐ 3-Function syringe LUZZANI
- ☐ 6-Function syringe LUZZANI
- ☐ Rack for micro blaster
 - ☐ with air connection ☐ Microetcher
 - ☐ without air connection ☐ Dento-Prep



Instruments positioned on chair base

- ☐ Yes
- ☐ No

Instruments


- ☐ Micromotor MX2 LED
- ☐ Turbine device 4VLM LUX
- ☐ Air pistol FREUDING

Quick-action coupling for external equipment

- ☐ Air/water
 - ☐ Water from bottle system
- ☐ Air

Instrument operation

- ☐ 4-Function pedal DYN
- ☐ 5-Function pedal Pedal
 - ☐ built in (standard)
 - ☐ free standing (optional)





3-Function pedal for operating special functions

- ☐ Yes
- ☐ No

Operating functions

- ☐ ORTHOLUX 200
- ☐ OP lamp ON/OFF
- ☐ special functions

ORTHODESK



Media supply ORTHODESK

- ☐ Yes (only possible, if no instrument is mounted on chair base)
☐ No

Media supply (2 instruments)

- ☐ Micromotor MX & MX2
☐ Micromotor MC3 IR
☐ Micromotor MC3 LED
☐ Turbine device 4VLM LUX
☐ 3-Function syringe LUZZANI
☐ 6-Function syringe LUZZANI

max. 1 instr.

Date:

Signature: